


The dispersing and deagglomeration of solids into liquids is an important application of ultrasonic devices. Ultrasonic cavitation generates high shear that breaks particle agglomerates into single dispersed particles.

The mixing of **powders into liquids** is a common step in the formulation of various products, such as [paint](#), [ink](#), [shampoo](#), [beverages](#), or [polishing media](#). The individual particles are held together by attraction forces of various physical and chemical nature, including van der Waals forces and liquid surface tension. This effect is stronger for higher viscosity liquids, such as polymers or resins. The attraction forces must be overcome in order to deagglomerate and disperse the particles into liquid media.

The application of mechanical stress breaks the particle agglomerates apart. Also, liquid is pressed between the particles. Different technologies are commonly used for the dispersing of powders into liquids. This includes high pressure homogenizers, agitator bead mills, impinging jet mills and rotor-stator-mixers.

High intensity ultrasonication is an interesting alternative to these technologies. When sonicating liquids the sound waves that propagate into the liquid media result in alternating high-pressure (compression) and low-pressure (rarefaction) cycles. This applies mechanical stress on the attracting electrostatic forces (e.g. van der Waals forces). **Ultrasonic cavitation** in liquids causes high speed liquid jets of up to 1000km/h (approx. 600mph). Such jets press liquid at high pressure between the particles and separate them from each other. Smaller particles are accelerated with the liquid jets and collide at high speeds. This makes ultrasound an **effective means for the dispersing** and deagglomeration but also for the [milling and fine grinding](#) of micron-size and [sub micron-size](#) particles.



www.imlab.com

Dispersing and Deagglomeration in Any Scale

Hielscher offers ultrasonic devices for the dispersing and deagglomeration of **any volume** for **batch or inline** processing. [Ultrasonic laboratory devices](#) are used for volumes from 1.5mL to approx. 2L. [Industrial ultrasonic devices](#) are used in the process development and production for batches from 0.5 to approx 2000L or flow rates from 0.1L to 20m³ per hour.

In the table below, you will find general device recommendations depending on the batch volume or flow rate to be processed. Click at the device type to get more information on each device.

Batch Volume	Flow Rate	Recommended Devices
0.5 to 1.5mL	n.a.	Vial Tweeter
1 to 500mL	10 to 200mL/min	UP100H
10 to 2000mL	20 to 400mL/min	UP200S , UP400S
0.1 to 20L	0.2 to 4L/min	UIP1000hd , UIP2000
10 to 100L	2 to 10L/min	UIP4000
n.a.	10 to 100L/min	UIP16000
n.a.	larger	cluster of UIP16000

Easy to Scale Up

Different from other dispersing technologies, ultrasonication can be **scaled up easily** from lab to production size. Laboratory tests will allow to select the required equipment size accurately. When used in final scale, the **process results are identical** to the lab results.

Robust and Easy to Clean



Ultrasonic power is transmitted into the liquid by the **sonotrode**. This is a typically rotary symmetric part, that is machined from solid Titanium. This is also the only moving/vibrating wetted part. It is the only part, that is subject to wear and it can be easily replaced within minutes. Oscillation-decoupling flanges allow to mount the sonotrode into open or closed pressurizable containers or flow cells in any orientation. No bearings are needed. All other wetted parts are generally made of stainless steel. Flow cell reactors have simple geometries and can **easily be disassembled** and wiped out. There are no small orifices or hidden corners.

Ultrasonic Cleaner in Place

www.imlab.com

Ultrasound is well known for its cleaning applications, such a surface, part cleaning. The ultrasonic intensity used for dispersing applications is much higher than for typical ultrasonic cleaning. When it comes to the cleaning of the wetted parts of the ultrasonic device, the ultrasonic power can be used to **assist cleaning** during flushing and rinsing, as the ultrasonic **cavitation removes particles** and liquid residues from the sonotrode and from the flow cell walls.