



Sieves cleaning –Nettoyage de tamis

Sieves cleaning module SRH 4/200

with integrated turning device for Elma ultrasonic units

- time-saving and cost-reducing intensive cleaning of analytical screens
- cleaning of up to 4 screens in one cleaning cycle (diameter of screens up to 200/203 mm (8"))
- analytical sieves of different makes can be cleaned at the same time

Only for baths S300H and S450H (one module per bath) and for bath S900 H (two modules per bath)

Référence	Libellé	P.U. HTVA
100 9133	Module avec rotation pour tamis, 230V	631.00



IMLAB sarl Centre d'Affaires de l'Horlogerie 48 rue des Canoniers 59000 LILLE	TEL : 03 20 55 19 11 FAX: 03 20 55 20 85	imlab@wanadoo.fr www.imlab.com
--	---	--