

**IMLAB sarl**  
**Instrumentation**  
**Matériels de laboratoire**

Centre d'Affaires de l'Horlogerie  
48 RUE DES CANONNIERS  
59000 LILLE



**03 20 55 19 11**

**03 20 55 20 85**

[imlab@wanadoo.fr](mailto:imlab@wanadoo.fr)

[www.imlab.com](http://www.imlab.com)

**Tarif HIELSCHER - 2009**

<b>Reference</b>	<b>Description</b>	<b>Price EUR HT</b>
<b>UP50H</b>	Ultrasonic processor for manual operation or stand use, 50 wats	2.200,00
<b>UP100H</b>	Ultrasonic processor for manual operation or stand use, 100 wats	2.534,00
<b>Sonotrode for UP50H and UP100H</b>		
<b>MS0.5</b>	Sonotrode made of titanium, diam 0,5 mm, lenght 80 mm, for samples from 0,01 ml to 0,5 ml	311,00
<b>MS1</b>	Sonotrode made of titanium, diam 1 mm, lenght 80 mm, for samples from 0,1 ml to 5 ml	311,00
<b>MS2</b>	Sonotrode made of titanium, diam 2 mm, lenght 80 mm, for samples from 2 ml to 50 ml	253,00
<b>MS3</b>	Sonotrode made of titanium, diam 3 mm, lenght 80 mm, for samples from 5 ml to 100 ml	253,00
<b>MS7</b>	Sonotrode made of titanium, diam 7 mm, lenght 80 mm, for samples from 10 ml to 250 ml	207,00
<b>Sonotrode for UP100H</b>		
<b>MS10</b>	Sonotrode made of titanium, diam 10 mm, lenght 80 mm, for samples from 20 ml to 500 ml	230,00
<b>Accessories for UP50H and UP100H</b>		
<b>ST1-16</b>	Stand, diam. 16mm height 600 mm, base 300 x 150mm	196,00
<b>SB2-16</b>	sound-protection box, with vertically adjustable table and diam 16 mm stand bar for UP 50H or UP100H	1.037,00

**Tarif HIELSCHER - 2009**

Reference	Description	Price EUR HT
<b>UP200S</b>	Ultrasonic processor for stand use, 200 wats	2.880,00
<b>UP200H</b>	Ultrasonic processor for manual operation or stand use, 200 wats	3.917,00
<b>Sonotrode for UP200S and UP200H</b>		
<b>S1</b>	Sonotrode made of titanium, diam 1 mm, lenght 100 mm, for samples from 0,1 ml to 5 ml	311,00
<b>S2</b>	Sonotrode made of titanium, diam 2 mm, lenght 100 mm, for samples from 2 ml to 50 ml	253,00
<b>S3</b>	Sonotrode made of titanium, diam 3 mm, lenght 100 mm, for samples from 5 ml to 200 ml	253,00
<b>S7</b>	Sonotrode made of titanium, diam 7 mm, lenght 100 mm, for samples from 20 ml to 500 ml	230,00
<b>S14</b>	Sonotrode made of titanium, diam 14 mm, lenght 100 mm, for samples from 50 ml to 1000 ml	230,00
<b>S40</b>	Sonotrode made of titanium, diam 40 mm, lenght 100 mm, for samples from 100 ml to 2000 ml	518,00
<b>Accessories</b>		
<b>ST1-16</b>	Stand, diam. 16mm height 600 mm, base 300 x 150mm	196,00
<b>SB1-16</b>	sound-protection box, with vertically adjustable table and diam 16 mm stand bar for UP 200S or UP400S	1.037,00
<b>SB3-16</b>	sound-protection box, with vertically adjustable table and diam 16 mm stand bar for UP 200H	1.037,00
<b>UP400S</b>	Ultrasonic processor for stand use, 400 wats	4.435,00
<b>Sonotrode for UP400S</b>		
<b>H3</b>	Sonotrode made of titanium, diam 3 mm, lenght 100 mm, for samples from 5 ml to 200 ml	294,00
<b>H7</b>	Sonotrode made of titanium, diam 7 mm, lenght 100 mm, for samples from 20 ml to 500 ml	294,00
<b>H14</b>	Sonotrode made of titanium, diam 14 mm, lenght 100 mm, for samples from 50 ml to 1000 ml	294,00
<b>H22</b>	Sonotrode made of titanium, diam 22 mm, lenght 100 mm, for samples from 100 ml to 2000 ml	294,00
<b>H40</b>	Sonotrode made of titanium, diam 40 mm, lenght 100 mm, for samples from 200 ml to 4000 ml	518,00
<b>Accessories</b>		
<b>ST1-16</b>	Stand, diam. 16mm height 600 mm, base 300 x 150mm	196,00
<b>SB1-16</b>	sound-protection box, with vertically adjustable table and diam 16 mm stand bar for UP 200S or UP400S	1.037,00

**Conditions de vente 2009**

Colonne Prix euro HT = prix de vente HORS TAXES, TVA 19,6% en sus

Livraison: franco si commande &gt; 400,00 EUR

Frais transport si commande &lt; 400,00 EUR: 25,00 EUR

Paiement: 30 jours date de facture

**GARANTIE: UNE ANNEE CONTRE TOUT VICE DE FABRICATION**