

Introduction Package: C-MAG HS 7



C-MAG HS 7 Package

C-MAG HS 7
+ ETS-D5

List price
~~749,00~~ EUR

Package price
549,00 EUR

Ident. No. 9015900

C-MAG HS 7
Magnetic stirrer with heating
new: with ceramic plate

- **Fixed safety circuit**
Shuts off heater if temperature rises above 550 °C
- **Hot top indicator**
Light flashes to warn user that the surface is hot to the touch
- **Digital display**
Allows for accurate and repeatable temperature settings
- **Powerful motor**
Stirs viscous materials with a maximum volume of 15 l (H₂O)
- **Strong magnet**
High magnetic adhesion prevents stir bar decoupling
- **Microprocessor control**
Sends feedback to the motor and hotplate to maintain constant speed and set temperature
- **Ceramic top**
Seamless one-piece glass-ceramic top provides excellent resistance to chemicals
- **Temperature controller**
Connection for contact thermometer enables precise temperature control of the medium
- **Elevated control panel**
Minimizes contact with spills and allows for easy reading and adjustment of controls



- + ETS-D5 Electronic contact thermometer
- + H 16 V Support rod
- + H 38 Holding rod
- + H 44 CLAMP

